


Disposable Surgical Pencil



Features

- Re-positionable stickers to ensure traceability in the patient's file
- Double sterilisation pouches peelable with sterilisation spot to answer to new regulation requirements ISO 11607-1
- Extra flexible cable to facilitate hand movements
- Hexagonal anti-rotation system blade for all electrodes 2.34 mm
- Compatibility of connectors VALLEY LAB ou ERBE type to guarantee risk free connections
- Certified quality - Class IIB

Description	Part Number	Connector	Cable length
Hand control - Blade electrode	BE 100	Valley Lab	3 m
Foot control - Blade electrode	BE 500	Valley Lab	3 m
Hand control - Blade electrode	SE 100	Valley Lab	3 m
Foot control - Blade electrode	SE 100 R.B	Valley Lab	3 m



JK Medical Systems Pvt.Ltd.

1068, Munusamy Salai, West K.K. Nagar, Chennai-600078, India
 Tel : 044 23664761 / 23664762 / 23664763 Fax : 044 23664764
 Email : sales@jkmedical.co.in Website : www.jkmedical.co.in



88 Avenue du General de Gaulle - 93541
 Bagnolet Cedex email : info@comepa.com
 website : www.comepa.com