

Bacterial Viral Filter

For Ventilator Use



Specification

- Description: Electrostatic Filter
- Dimension: H - 67.2mm, W - 68.5mm
- Filtration efficiency: BFE 99.9999%; VFE 99.9999%
- Connection inlet: Patient side 22M/15F machine side 22F/15M
- Weight: 23g
- Applications: Anaesthesia, ICU, Emergency, Homecare. The filter can be placed patient side or machine side

Order Information

Cat No	Description	Pkts
4222/701	BV Filter	50 Nos.

CE
0473

Marketed & Distributed by



JK Medical Systems Private Limited

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in



GVS Filter Technology UK Ltd.

NFC House, Mellishaw Lane | Morecambe | Lancashire | UK LA3 3EN | www.gvs.com