



CE 0123

EOLO is a piston compressor nebulizer for aerosoltherapy. Designed for home-care with convenient handle and a large compartment for easy transport of unit, accessories and medications.

### STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
- Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)
- Air Tube
- Air Filter

### TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump
Power Feeding	230V-50Hz - other configurations on request
Power consumption	170 VA
Max Pressure	2.5 bar      36 psi      250 kPa
Operating Pressure	1.1 bar      16 psi      110 kPa
Max Air Flow	14 l/min
Operating Air Flow	5 l/min
Neb-rate with 4 ml solution	0,40 ml/min
Mass Median Aerodynamic Diameter	2,44 $\mu$ m (MMAD) according to EN-13544-1
Respirable Fraction (< 5 $\mu$ m)	> 80%
Duty Cycle	Non-stop operation
Noise Level (EN 13544-1 Standard)	Approx. 55 dBA
Weight	1,65 kg
Size (Main Unit)	cm 14 x 12 x 22
Years of Warranty	1
Shipping carton	6
Place of Manufacturing	Italy

Italian  
Medical  
Touch  
MADE IN ITALY

CM005.REV1 - EDITION NOVEMBER 2014

### CODE

RE 300400 ————— EOLO