



CE 0123

MIKO is a piston compressor nebulizer for aerosoltherapy. Designed for home-care with compact size and carrying bag for easy transport of unit, accessories and medications. Recently upgraded for non-stop operation in case intensive treatments are required.

STANDARD CONSUMABLES INCLUDED

- HI-FLO Jet Nebulizer
- Adult Masks with elastic strings
- Child Masks with elastic strings
- Mouthpiece
- Nosepiece (non invasive)
- Air Tube
- Air Filter
- Carrying Bag (MIKO model only)

TECHNICAL FEATURES

Compressor	Oilless and maintenance-free piston pump		
Power Feeding	230V-50Hz - other configurations on request		
Power consumption	170 VA		
Max Pressure	2.2 bar	32 psi	220 kPa
Operating Pressure	0.95 bar	13.8 psi	95 kPa
Max Air Flow	12 l/min		
Operating Air Flow	4,8 l/min		
Neb-rate with 4 ml solution	0,35 ml/min		
Mass Median Aerodynamic Diameter	2,44 μ m (MMAD) according to EN-13544-1		
Respirable Fraction (< 5 μ m)	> 80%		
Duty Cycle	Non-stop operation		
Noise Level (EN 13544-1 Standard)	Approx. 56 dBA		
Weight	1,45 kg		
Size (Main Unit)	cm 14 x 17 x 11		
Years of Warranty	1		
Shipping carton	6		
Place of Manufacturing	Italy		

Italian
Medical
Touch
MADE IN ITALY

CM005.REV1 - EDITION NOVEMBER 2014

CODE

- RE 300600/03 ————— MIKO- with carrying bag
 RE 300600/12 ————— MIKO BASIC - without carrying bag



JK Medical Systems Pvt.Ltd.

1068, Munusamy Salai, West K.K. Nagar, Chennai-600078, India
 Tel : 044 23664761 / 23664762 / 23664763 Fax : 044 23664764
 Email : sales@jkmedical.co.in Website : www.jkmedical.co.in



via Ugo La Malfa 13 - 43010 Pilastrò (PR) - Italy
 Tel. +39 0521 637133/631138 Fax +39 0521 639041
 export@ca-mi.it - vendite@ca-mi.it - www.ca-mi.it