

PARI COMPACT

Quality Compressor Nebulizer



- Inhalation device for diseases of the lower airways, such as COPD, Asthma or Bronchitis.
- Proven efficiency: Includes PARI LC PLUS nebulizer with mouthpiece and child's mask.

Technical data

Operation	1.0 bar
Weight	1.2 Kg
Dimension (wxhxd)	10x14x16(cm)

Aerosol Performances

TOR [mg/min]	418
MMD [μm]	3.9
Mass percentage $<5\mu\text{m}$	64%

Measured with the Malven MaserSize X (according to Fraunhofer-Model) at 21° C and 50% relative humidity. Nebulized medium: 0.9% NaCl (5ml). Inspiratory flow 20 l/min

PARI COMPACT and its Accessories

 <p>PARI COMPACT Inhalation device 230-240 V ~ /50Hz with PARI LC PLUS® Nebulizer, Child mask, Mouthpiece universal and Connecting tubing Item No. 052G1002</p>	 <p>PARI LC PLUS® Consists of LC PLUS® Nebulizer, Connection Tubing (1.2m) and Filter Item No. 022G8010</p>	 <p>Filter Holder with Filter Item No. 041B6500</p>
 <p>PARI Child Mask Soft "Spiggy" with Elastic Band Item No. 041G0741</p>	 <p>PARI Adult Mask Soft with Elastic Band Item No. 041G0740</p>	 <p>Connection Tubing (1.2 m, f/f) Item No. 041B4592</p>

*Specifications subject to change without prior notice.



JK Medical Systems Pvt Ltd

1068, Munusamy Salai, West K K Nagar, Chennai - 600078, India
Tel : 044 23664761/23664762 | sales@jkmedical.co.in | www.jkmedical.co.in



PARI GmbH

Moosstrasse 3 D-82319 Starnberg
Tel : +49 (0) 8151 279-0 | Fax: +49 (0) 8151 279-101
www.pari.de