

# PARI VELOX

## Nebuliser for mobile use



Velox



Velox Junior

### Description

Mobile nebuliser device based on vibrating membrane technology. Suitable for treatment of acute or chronic diseases of the lower airways. The ergonomic design and small number of components makes it easy to use, and the extremely quiet operation guarantees discreet, relaxed treatment. With optional mains or battery operation, small size and light weight, the VELOX is also an ideal device for use while on travel.

### Technical data

Inhalation time	2–4 minutes (for 2.5 ml)
Power sources	Mains power: Multi-voltage power adapter (100–240V~/50-60Hz) for use all over the world Battery power: 3 × 1.5 V Mignon AA LR6/alkaline batteries or 3 × 1.2 V NiMH rechargeable batteries Note: The device does not have a battery charging function. Operating time on battery power: With new, high-power batteries from a trusted brand, the operating time is approx. 4 hours. New rechargeable batteries (AA), each with 2000 mAh energy storage capacity, achieve comparable operating times.
Weight	110g (without batteries) 140g (with batteries)

### Aerosol Characteristics

Total Output Rate	497 mg/min
MMD	3.8 µm
Mass fractions < 5 µm	74 %



**JK Medical Systems Pvt Ltd**

1068, Munusamy Salai, West K K Nagar, Chennai - 600078, India  
Tel : 044 23664761/23664762 | sales@jkmedical.co.in | www.jkmedical.co.in



**PARI GmbH**

Moosstrasse 3 D-82319 Starnberg  
Tel : +49 (0) 8151 279-0 | Fax: +49 (0) 8151 279-101  
www.pari.de