

PARI BOY[®] mobile S Nebulizer

For True Portability



Features

- The portable compressor nebuliser with advanced rechargeable battery - for treatment anytime, anywhere.
- The lightweight member of the PARI family of nebuliser systems-small, powerful compressor provides efficient therapy.
- When battery is fully charged, the device can be used upto 50 minutes, equal to about 8 inhalation sessions. (with a medication fill volume of 2.5ml).
- Supplied with a 12 volt adapter and universal mains charging adapter (100 to 240 V / 50-60Hz), the PARI BOY mobile S can be used anywhere in the world.
- For use in respiratory tract disorder such as asthma, bronchitis, chronic obstructive pulmonary diseases, or cystic fibrosis.
- Your constant companion - for work, leisure, or trips abroad.
- The battery can be completely recharged in about 2 hours.
- Suitable for use by everyone aged 4 and over.

Standard Supplies

PARI BOY mobile S compressor with rechargeable battery, 12 volt adapter, mains charging adapter (100 to 240 V/50-60 Hz, PARI LC SPRINT nebulizer, connecting tube, mouthpiece, mask and carrying case.

Aerosol Characteristics

MMD [μm]	3.8*
TOR [mg/min]	440
%<5 μm	65

* Inspiratory flow 20 l/min, measured by Malvern MasterSizer X at 23 °C and 50% relative humidity. Nebulised medium: 0.9% NaCl (5 ml).

Technical data for the compressor

Pressure	1.0 bar
Electrical connection	100 - 240 V/50-60Hz
Dimensions (W x H x D) [cm]	8.9 x 4.5 x 14.8
Battery	Lithium ion, rechargeable
Weight [g]	520
Warranty	1 year

*Specifications subject to change without prior notice.



JK Medical Systems Pvt Ltd

1068, Munusamy Salai, West K K Nagar, Chennai - 600078, India
Tel : 044 23664761/23664762 | sales@jkmedical.co.in | www.jkmedical.co.in



PARI GmbH

Moosstrasse 3 D-82319 Starnberg
Tel : +49 (0) 8151 279-0 | Fax: +49 (0) 8151 279-101
www.pari.de