

SpiroSafe Filter

For Spirometry Use



Description

Disposable bacterial/viral filter for pulmonary function test comes with 3 different fittings at the machine end and practically fits into all machines without a need of an adaptor. Ergonomically designed with high efficiency filtration of 99.99% of bacteria and virus.

Specification

- Housing: White polypropylene
- Electrostatic Spirometry Filter
- 99.99% bacterial/vital protection
- Weight: 10.87gms
- Height: 8.5 cm
- Resistance 0.7cm H 20/L/s @720 L/min
- Dead volume 55ml
- Latex Free

Order Information

Cat No	SpiroSafe Small	SpiroSafe Medium	SpiroSafe Large	Pkts
FF010010	●			50 Nos.
FF010011		●		50 Nos.
FF010012			●	50 Nos.



JK Medical Systems Private Limited

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in