


Disposable Surgical Pencil



Features

- Re-positionable stickers to ensure traceability in the patient's file
- Double sterilisation pouches peelable with sterilisation spot to answer to new regulation requirements ISO 11607-1
- Extra flexible cable to facilitate hand movements
- Hexagonal anti-rotation system blade for all electrodes 2.34 mm
- Compatibility of connectors VALLEY LAB ou ERBE type to guarantee risk free connections
- Certified quality - Class IIB

Description	Part Number	Connector	Cable length
Hand control - Blade electrode	BE 100	Valley Lab	3 m
Foot control - Blade electrode	BE 500	Valley Lab	3 m
Hand control - Blade electrode	SE 100	Valley Lab	3 m
Foot control - Blade electrode	SE 100 R.B	Valley Lab	3 m

Distributed and Marketed By:



JK Medical Systems Private Limited

1068, Munusamy Salai, West K K Nagar, Chennai-600078, India

Tel : 044 23664761 / 23664762 | sales@jkmedical.co.in | www.jkmedical.co.in



88 Avenue du General de Gaulle - 93541

Bagnolet Cedex email : info@comepa.com

website : www.comepa.com