

Disposable Electrosurgical Plate



Features

- Traceability** : The re-positionable sticker ensures the direct traceability in the patient's file.
- Adaptability** : The large range of shapes and sizes available permits to choose the best product following the positioning and morphology of the patient.
- Safety** : The low impedance of specially developed hydrogel ensures a rise of the skin temperature less than 6° (it is the limit accepted by AAMI HF-18, 2001 regulation).
- Components** : Conductive hydrogel with high adhesive and hydrophilic power. PE Foam support or non woven support impervious to liquid infiltration - Conductive medical aluminium.
- Compatibility** : All PLACOM plates have been certified compatible with electrosurgery units and comply with European regulations. They are certified CE0459, Class IIB

Description	Part Number	Connector
Prewired 2-Part Plates (Adult)	MI 3.05.2500 C	Valley Lab
Prewired 2-Part Plates (Child)	MI 3.05.2600 C	Valley Lab
2-Part Plates (Adult)	MF 3.05.3000.01	Any Machine
2-Part Plates (Child)	MF 3.05.3000.02	Any Machine
1-Part Plates (Adult)	MF 3.05.2000.01	Any Machine
1-Part Plates (Child)	MF 3.05.2000.01	Any Machine

Distributed and Marketed By:



JK Medical Systems Private Limited

1068, Munusamy Salai, West K K Nagar, Chennai-600078, India

Tel : 044 23664761 / 23664762 | sales@jkmedical.co.in | www.jkmedical.co.in



88 Avenue du General de Gaulle - 93541

Bagnolet Cedex email : info@comepa.com

website : www.comepa.com