


IP - 1010

Tabletop pulse oximeter



Features

- Easy to transport design
- Nellcor compatible SpO₂
- Protection from motion artifact
- Accurate during a low perfusion
- SpO₂ averaging (4/8/12/16) technology
- Up to 8 hours use with an inbuilt rechargeable battery / Light weight 1.3 kg
- Wide display with Selective Plethmographic mode, SpO₂ Averaging, Trend Data and Numeric Mode

Standard Specifications

SpO₂

SpO ₂ Range	0 ~ 100%
SpO ₂ Accuracy	70 to 100% $\pm 2\%$ / 50 to 69% $\pm 3\%$
SpO ₂ Averaging	4/8/12/16
Pulse rate Range	20 - 250 bpm
Pulse rate Accuracy	20 - 250 bpm $\pm 2\%$
Perfusion Range	0.05% - 20%

Alarm Limit

SpO ₂	High range 36% ~ 100% : Adjustment per 1% Low range 35% ~ 99% : Adjustment per 1%
Pulse rate	High range 21 ~ 250 BPM : Adjustment per 1 BPM Low range 20 ~ 249 BPM : Adjustment per 1 BPM
Alarm silence option	30/60/90/120 seconds duration

Power

Adaptor	Input 100 to 240V 50/60Hz, Output DC 12V, 1.5A
Power Consumption	30W
Battery Type	Li-Ion
Charging time	5 hours
Operating time	10 hours

Accessories

Standard	Adult SpO ₂ probe Optional: Paediatric and Neonatal probe
----------	---

Marketed & Distributed by:



JK Medical Systems Private Limited

1068, Munusamy Salai, West K.K. Nagar, Chennai-600078, India
Tel : 044 23664761 / 23664762 | sales@jkmedical.co.in
www.jkmedical.co.in



Infunix Technology Co., Ltd.

Seokjin B/D 2F, 11, Bongeunsa-ro 69-gil, Gangnam-gu, Seoul,
Korea, 135-868 | sales@infunix.com
www.infunix.com