


Spiroflo® II mouthpiece

For Spirometry Use



Description

Uniquely designed keeping in mind it's predecessor Spiroflo mouthpiece. Fits easily to patients of all ages without any leakage. Ergonomically designed for maintaining better flow rate.

Specification

- Height : 7cm
- Weight : 5.7gms
- Material : PE

Order Information

Cat No	Description	Min. Order Qty
SM010007	Spiroflo mouthpiece	100 nos
SM010004	Spiroflo mouthpiece - Jaeger m/c	100 nos



JK Medical Systems Private Limited

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in