

# TryFlo

## Incentive Spirometer



- Tryflo incentive spirometer is ideal for developing, improving and maintaining your respiratory fitness.
- Tryflo incentive is a deep breathing exerciser that has been scientifically constructed as a means of encouraging you to take a slow sustained maximal inspiration (SMI).
- Tryflo incentive spirometer is designed to help you perform normal breathing exercises. In turn, the spirometer will help you to prevent the possibility of respiratory complications or to improve your breathing pattern if you have a respiratory condition.

### Technical Specifications

Flow Rate	600 ml/sec (Red)/ 900 ml/sec (Yellow)/ 1200 ml/sec (Blue)
Dimension	135 x 135 x 68mm
Weight	150 gms

### Order Information

Cat No	Pkts
TF010070	24 Nos/ box



**JK Medical Systems Private Limited**

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in  
www.jkmedical.co.in