


# Surgical Instruments



## Bipolar Lap Instruments

- Bipolar laparoscopic instruments 5mm x 330mm for micro and macro MIS applications
- Different shape tips and 360° rotatable handle (Rotaris)
- Flat plug for ESU connectivity.



## Bipolar Scissors & Clamps

- Bipolar scissors - Single pin connectivity, one blade with ceramic insert, very thick insulation.
- Bipolar clamps - Coagulation effected by flat clamp / manual cutting by scissors incorporated, very thick insulation.



**JK Medical Systems Private Limited**

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India

Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in



**MICROMED Medizintechnik GmbH**

Eisenbahnstrabe 84

78573 Wurmlingen, Germany