

Surgical Instruments



Leyla Retractors

- A very dependable brain retractor system made of SS of highest hardness to withstand wear and tear
- Flexible cable and locking system for ease of usage
- Flexible Nitinol spatulas that retain original shape when heated



Vessel Clips

- Aneurysm clips available as standard and mini, permanent and temporary
- Different sizes, shapes and suitable applicators
- Applicators with shaft made of heat sensitive Nitinol
- Clips available in titanium and phynox
- MRI compatible upto 6 tesla



Container Systems

- High quality aluminium containers for storage of various instruments and instrument sets
- Option of perforation at bottom or top or both
- Different types of filters for single or multiple use
- Silicon beading assures sixty days sterility
- Silicon mats for micro instruments



JK Medical Systems Private Limited

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India

Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in

LAWTON

LAWTON GmbH & Co. KG

Wurtemberger Str. 23

78567 Fridingen, Germany