

Surgical Instruments



Plastic Surgery

- Very high precision micro surgery instruments range
- Full range available for plastic, micro vascular, hand reconstructive surgeries
- Vascular clips & approximators, arteriotomy forceps, nerve and tendon forceps offered individually or as sets with sterilization trays.



Vessel Sealer

- Specially designed Bipolar laparoscopy instrument, can be used for sealing blood vessels upto 7mm and modify bundle of tissue
- An in-built mechanism with a detachable blade to cut, avoiding change of instruments
- Luer port for easy cleaning of the shaft
- Comes with flat type Bipolar connectivity.



JK Medical Systems Private Limited

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India

Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in



MICROMED Medizintechnik GmbH

Eisenbahnstrabe 84

78573 Wurmlingen, Germany