

# Surgical Instruments




## Monopolar Tips

- Disposable and Reusable electrode tips for diathermy pencils
- Blade, knife, needle, ball, cup, loop and ribbon shaped tips, angled and straight.



## Bipolar Forceps

- High quality bipolar forceps with long lasting insulation
- Available in various tip sizes
- Straight, angled, bayonet shaped forceps
- Available as regular, non-stick, Argentum and with TC inserts



## ESU Cables

- Autoclavable bipolar and mono polar connecting cables to suit any ESU
- Copper strands coated with silver and double insulated
- Ensures perfect connectivity
- Available for both monopolar and bipolar hand instruments



**JK Medical Systems Private Limited**

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India

Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in



**MICROMED Medizintechnik GmbH**

Eisenbahnstrabe 84

78573 Wurmlingen, Germany