



JK Medical Systems Pvt Ltd
www.jkmedical.co.in

HME Filter (Paed)

For Ventilator Use



Specification

- Description: Electrostatic Filter with Heat Moisture Exchanger (HMEF)
- Dimension: H - 73.0mm, W - 48.0mm
- Filtration efficiency: BFE 99.9999%; VFE 99.9999%
- Connection inlet: Connections (ISO 5356-1)
- Weight: 19g
- Applications: Anaesthesia, ICU, Emergency, Homecare. The filter should be placed patient side between the breathing circuit and the catheter mount

Order Information

Cat No	Description	Pkts
9066/711	HME Filter	50 Nos.



Marketed & Distributed by



JK Medical Systems Private Limited

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in



GVS Filter Technology UK Ltd.

NFC House, Mellishaw Lane | Morecambe | Lancashire | UK LA3 3EN | www.gvs.com