

# HME Filter

For Ventilator Use



## Specification

- Description: Heat Moisture Exchanger Filter (HMEF) Clinic Clean blister packed
- Dimension: H - 77.0mm, W - 68.5mm
- Filtration efficiency: BFE 99.9999%; VFE 99.9999%
- Connection inlet: Connections (ISO 5356-1)
- Weight: 27g
- Applications: Anaesthesia, ICU, Emergency, Homecare. The filter should be placed patient side between the breathing circuit and the catheter mount

## Order Information

Cat No	Description	Pkts
4333/711	HME Filter	50 Nos.

CE  
0473

Marketed & Distributed by



**JK Medical Systems Private Limited**

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in



**GVS Filter Technology UK Ltd.**

NFC House, Mellishaw Lane | Morecambe | Lancashire | UK LA3 3EN | www.gvs.com