

LungChek

Peak Flow Meter



- Large, high visibility scale
- Available with European/ATS scale
- Movable asthma management color zones
- Customization possible for promotional purposes
- Meets and exceeds all international performance standards
- Available Colors: Purple, Blue & Green

Technical Specifications

Range	60-800 l/min
Scale Increment	10 l/min
Accuracy	<+/- 5% or 10 l/min
Repeatability	<+/- 3%
Calibration	100% (individual unit)
Weight	60 gms
Dimension	Height - 150mm Diameter - 42mm

Order Information

Cat No	Pkts
LC010050	01 No.



JK Medical Systems Private Limited

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in