

PeakFlo

mouthpiece



Description

- Made out of medical grade LLDPE material
- Adult & Paed use
- Available with one way valve (optional)

Specification

- Height : 62.5 mm
- Weight : 6 gms

Order Information

Cat No	Pkts
PF010077	50 Nos.



JK Medical Systems Private Limited

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in