

Spiroflo® Filter

For Spirometry Use



Description

Disposable filters for pulmonary function test with electrostatically charged filter medium for high efficiency filtration of 99.999% of virus and bacteria. Can be incorporated with patient friendly mouthpiece with negligible resistance for accurate measurements. Available in combo packs with mouthpiece and nose clip.

Specification

- Housing: White polypropylene
- Electrostatic Spirometry Filter
- 99.99% bacterial/viral protection
- Weight: 27gms
- Height: 8.5cm
- Resistance 0.7cm H 20/L/s @720 L/min
- Dead volume 75ml
- Latex Free

Order Information

Cat No	Description	Pkts
SF010012	Spiroflo Filter	100 Nos.
SF010014	Spiroflo Filter with mouthpiece	100 Nos
SF010015	Spiroflo Filter with mouthpiece & Nose Clip	100 Nos



JK Medical Systems Private Limited

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in