

Spiroflo® II mouthpiece

For Spirometry Use



Description

Uniquely designed keeping in mind its predecessor Spiroflo mouthpiece. Fits easily to patients of all ages without any leakage. Ergonomically designed for maintaining better flow rate.

Specification

- Height : 7cm
- Weight : 5.7gms
- Material : PE

Order Information

Cat No	Pkts
SM010007	50 Nos.



JK Medical Systems Private Limited

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in