


# Spiroguard Filter

For Spirometry Use



## Specification

- Description: Electrostatic Filter Clinic Clean pouch packed
- Dimension: H - 92.65 mm; W - 96.7 mm
- Filtration efficiency: BFE 99.9999%; VFE 99.9999% up to 0,027  $\mu\text{m}$
- Materials of construction: Filter media - Electrostatic; Housing - Polypropylene
- Weight: 36g
- Applications: lung functional test for use in Respiratory Pathophysiology, ICU and Rehabilitation, Sports Medicine, Occupational Medicine

## Order Information

Cat No	Description	Pkts
2800/21	Spiroguard Filter	50 Nos.

Marketed & Distributed by



**JK Medical Systems Private Limited**

(AN ISO 13485:2012 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in



**GVS Filter Technology UK Ltd.**

NFC House, Mellishaw Lane | Morecambe | Lancashire | UK LA3 3EN | www.gvs.com