


Easyflo Bacterial Viral Filter

For Ventilator Use



Specification

- Description: Electrostatic Filter
- Dimension: H - 63mm, W - 68mm
- Filtration efficiency: BFE 99.99990%; VFE 99.99975%
- Connection inlet: Patient side 22M/15F machine side 22F/15M
- Weight: 21g
- Applications: Anaesthesia, ICU, Emergency, Homecare. The filter can be placed patient side or machine side

Order Information

Cat No	Description	Qty
7010	BV Filter - Adult	100 Nos
7011	BV Filter - Paed	100 Nos



JK Medical Systems Private Limited

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in