

Filteratte Pro

For ndd spirometry Use



Description

A specially designed bacterial/viral filter for use with NDD PFT machines only. It fits universally for adult & paed use. Comes with high efficiency filtration of BEF is 99.9992% and VEF is 99.9970%.

Specification

- Housing: White Polyethylene
- Electrostatic Spirometry Filter
- Dimensions: H- 149.5mm, W- 56.5mm
- BEF is 99.9992% & VEF is 99.9970%
- Resistance < 1.6cmH₂O @ 32LPM
- Dead volume 50ml
- Weight: 16.17 gms

Order Information

Cat No	Description	MOQ
5052	Filteratte Pro	200 nos



JK Medical Systems Private Limited

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in

www.jkmedical.co.in