


Easyflo HME Filter

For Ventilator Use



Specification

- Description: Heat Moisture Exchanger Filter (HMEF) with sterile pack
- Dimension: H - 73mm, W - 68mm
- Filtration efficiency: BFE 99.9999%; VFE 99.9975%
- Connection inlet: Connections (ISO 5356-1)
- Weight: 23g
- Applications: Anaesthesia, ICU, Emergency, Homecare. The filter should be placed patient side between the breathing circuit and the catheter mount

Order Information

Cat No	Description	Qty
7020	HME Filter - Adult	100 Nos
7021	HME Filter - Paed	100 Nos



JK Medical Systems Private Limited

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in