

Consumables



Indspire Nebulizer Kit

Comes in a set of mask (Adult/Paed), Chamber and Tubing. Also available with T-piece mouthpiece combination

Specification

Total Output Rate	450mg/min
MMD	4.2µm
<5µm%	63%

Order Information

Cat No	Kit with Adult Mask	Kit with Paed Mask	Kit with Mouthpiece	Kit with Adult & Paed Mask	QTY
2001	•				1 No.
2002		•			1 No.
2003			•		1 No.
2004				•	1 No.



Oxygen Mask

Comes in a set of mask (Adult/Paed), with 2 metre tubing.

Order Information

Cat No	Description	Pkts
8051	Adult mask Set	50 Nos
8052	Paed mask Set	50 Nos



High Concentration Mask

It is used for the application of heated and humidified blended air and oxygen at high flow rates. This high flow can provide more constant inspiratory oxygen concentrations than conventional oxygen therapy and can also generate some positive end expiratory pressure

Order Information

Cat No	Description	Pkts
8021	Adult mask Set	50 Nos
8022	Paed mask Set	50 Nos



JK Medical Systems Private Limited

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in