

SpiroSafe Filter

For Spirometry Use



Description

Disposable bacterial/viral filter for pulmonary function test comes with 3 different fittings at the machine end and practically fits into all machines without a need of an adaptor. Ergonomically designed with high efficiency filtration of BEF is 99.9992% and VEF is 99.9970%.

Specification

- Housing: White polypropylene
- Electrostatic Spirometry Filter
- BEF is 99.9992% & VEF is 99.9970%
- Weight: 10.87gms
- Height: 8.5 cm
- Resistance < 1.6cmH₂O @ 32LPM
- Dead volume 55ml
- Latex Free

Order Information

Cat No	Description	MOQ
5000/S	SpiroSafe Filter (Small) - NDD m/c	200 nos
5000/M	SpiroSafe Filter (Medium) - MIR & Cosmed m/c	200 nos
5000/L	SpiroSafe Filter (Large) - Jaeger m/c	200 nos



JK Medical Systems Private Limited

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in
www.jkmedical.co.in