

# TryFlo

## Incentive Spirometer



- Tryflo incentive spirometer is ideal for developing, improving and maintaining your respiratory fitness.
- Tryflo incentive is a deep breathing exerciser that has been scientifically constructed as a means of encouraging you to take a slow sustained maximal inspiration (SMI).
- Tryflo incentive spirometer is designed to help you perform normal breathing exercises. In turn, the spirometer will help you prevent the possibility of respiratory complications or to improve your breathing pattern if you have a respiratory condition.

### Technical Specifications

|           |   |
|-----------|---|
| Flow Rate | 600 ml/sec (Red)/ 900 ml/sec (Yellow)/ 1200 ml/sec (Blue) |
| Dimension | 135 x 135 x 68mm  |
| Weight    | 150 gms   |

### Order Information

| Cat No | Qty         |
|--------|-------------|
| 6070   | 10 Nos/ box |



**JK Medical Systems Private Limited**

(AN ISO 13485:2016 & WHO-GMP CERTIFIED COMPANY)

1068, Munusamy Salai, West KK Nagar, Chennai-600 078, India | Ph:044-23664761/ 62 | sales@jkmedical.co.in  
www.jkmedical.co.in